



- Since 1980
- Automatic Compact Blistering Machines
- R&D Blistering Machines
- Deblistering machines (Automatic and Semi-automatic)
- Format Parts for high Speed Blistering Machines
- Service and Maintenance

Milan - Italy

4 Description of the Machine

The Automatic Deblistering Machine model SIMPATY F160 presents itself to our customers with accurate design and once again evolved from the first version. The clean lines draw a machine that moves on a support with wheels suitable for handling in the production departments. Technology of the automatic loader placed on the head of the machine allows to accelerate the production up to a blister at second always relative to the type of the product and the length of the format. The recovery of the product is done with the positioning and the set up of the pressure rollers and blades to engrave and cut the covering material. Subsequently the punches are positioned in correspondence of each pocket and the product is dropped into containers of stainless steel for a safe recovery while the blisters are recovered in another container, both positioned within the design of the F160. Today, on the new version it is also possible to put dedicated format ready to be mounted and dedicated for each blister to avoid the positioning of each deblistering accessory to have faster set up.



wondershare™

PDF Editor

- Since 1980
- Automatic Compact Blistering Machines
- R&D Blistering Machines
- Deblistering machines (Automatic and Semi-automatic)
- Format Parts for high Speed Blistering Machines
- Service and Maintenance

5 Machine's Layout

MACHINE'S DIMENSIONS (LxPxH)

PACKAGING DIMENSIONS (LxPxH)

WEIGHT OF THE MACHINE (WITHOUT PACKAGING)

VOLTAGE

LIGHT INSIDE WORKING AREA

940x560xH1630MM

1170x670xH1810MM

120 KG

230V 50 Hz

NOT SUPPLIED

