



Handling



Dispensing



Sieving



Blending



Container



Cleaning

High Containment System

Flexible und precise



SERVOLIFT

Your specialist for
complete handling-, blending-
and cleaning technologies

SERVOLIFT GmbH
Albert-Einstein-Straße 9
77656 Offenburg Germany

T. +49 (0) 781 6100 0
F. +49 (0) 781 6100 400

info@servolift.de
www.servolift.de

► Flexible High Containment System

SERVOLIFT offers a solution for another challenging task frequently found in the field of high-containment applications: Milling of pharmaceutical powders by using only mobile handling equipment.

All heights are safely held by mechanical height stops - these together with the positioning frame guarantee repeatability of all positions.



Small containers and also drums can be used

ConiWitt milling-head

Retractable height stop

High-Containment split-valve system

Easy and smooth drum docking by manually lifting

Manually adjustable positioning frame

Features:

- Entire system complies with latest GMP regulations
- Mobile system allows a maximum of flexibility
- High-containment milling and product-transfer of small batches
- Milling-head can be easily removed
- Containers and drums per customers application
- Design according to customers specification