



BLISTER SOLIDA

Tablets  
Oblongs  
Capsules

## BS-600

## BLISTER MACHINE | SOLIDA

The BS-600 combines innovative design with high flexibility and achieves a reliable output up to 600 blisters per minute. Depending on the specific application the machine can be equipped with any combination of several options. The machine is available with rotary or platen sealing.

The versatile machine concept is designed for thermoforming and cold forming processes.



### TECHNICAL INFO

max. output	600 blisters/min.
max. forming area	258 x 205 mm, optional 258 x 290 mm
max. forming depth	12 mm, optional 20 mm



Dedicated high speed feeding system



Combi feeding system with channels