



I Holland Segment Cleaning

Easy, repeatable cleaning process



I Holland brings to you the fastest and most efficient method for cleaning segments. Our specially designed segment fixtures fit perfectly into the ultrasonic tank and leave the bores unobstructed, meaning effective cleaning of the entire segment. The pivot design allows for an easy drying and drainage process, and the basket coating prevents scratching and damage... You don't have to worry about a thing.

On top of all of this, our new ultrasonic machine cleans segments, punches and dies, saving you valuable space. The ultrasonic method provides repeatable cleaning results compared to manual cleaning so you can remove stubborn contaminants from all of your tooling in a lot less time!

Cleaning method	Time taken for a full set of tools (Fette 3090i)		Total time
	220 punches	5 segments	
I Holland Segment Cleaner	10'	50'	1h
Manual Cleaning	220'	150'	6h 10'
Time saved per set	1h 50'	1h 40'	5h 10'

www.tablettingscience.co.uk

I Holland Limited, Meadow Lane, Long Eaton, Nottingham. NG10 2GD. England
 t: +44 (0) 115 972 6153 // f: +44 (0) 115 973 1789 // e: info@iholland.co.uk